

Tyvek® Offers Unmatched Sterilization Compatibility

Article | February 19, 2018



Unlike medical-grade papers and films, DuPont™ Tyvek® offers sterilization compatibility with all of the most commonly used methods for sterilizing medical devices. These include: [ethylene oxide \(EO\)](#), [gamma](#), [electron-beam](#), [hydrogen peroxide gas plasma](#) and [steam](#) (under controlled conditions). That's because Tyvek® is made of 100% high-density polyethylene (HDPE), which is extremely stable when exposed to sterilant gases and high-energy sterilization processes. What's more, Tyvek® is specially engineered to enable sterilant gases and steam to penetrate and escape quickly.

When it comes to sterilization compatibility, Tyvek® offers the broadest range of options. And no matter which process you choose—EO, gamma, electron-beam, steam (under controlled conditions) or low-temperature oxidative sterilization processes—Tyvek® will retain its superior protective properties of microbial barrier and strength, as well as its color and flexibility.

Material compatibility with various sterilization methods

	DuPont™ Tyvek®	Coated, latex saturated medical-grade paper	Medical film
Ethylene oxide (EO)	Yes	Yes	No
Gamma radiation	Yes	Yes	Yes
Electron-beam radiation	Yes	Yes	Yes
Steam	Yes ¹	Yes ²	No
STERRAD®	Yes	No	No

1. Under controlled conditions (250°F to 260°F [121°C to 127°C]) at 30 psi for 30 minutes.
2. May become brittle.

Ethylene Oxide [EO]

[Ethylene oxide \(EO\)](#) does not readily adsorb on Tyvek® and is released more rapidly than from cellulosic materials such as medical-grade papers, including synthetic fiber-reinforced paper. Tyvek® maintains its superior strength and microbial barrier properties after EO sterilization.

Gamma Radiation

After exposure to [gamma radiation](#) up to 100 kGy, Tyvek® maintains its superior microbial barrier and the impact on strength properties is limited. These properties are also maintained after irradiation

followed by exposure to accelerated and real-time aging.

Electron-Beam Radiation

After exposure to **electron-beam radiation** up to 100 kGy, Tyvek® maintains its superior microbial barrier and strength properties. Although specific studies of electron-beam sterilization followed by aging have not been conducted, we have not seen any effects other than those shown after **gamma radiation** and aging. This is because HDPE is radiation stable.

Hydrogen Peroxide Gas Plasma

Tyvek® 4057B is suitable for use with the STERRAD® Sterilization System from Advanced Sterilization Products (ASP), Division of Ethicon Inc., a Johnson & Johnson company. This sterilization method uses **low-temperature gas plasma** to enable sterilization of heat-labile devices.

Caution should be used when choosing methodologies for package integrity testing after multiple doses of low-temperature oxidative sterilization because the water resistance of the material can be altered.

Medical-grade papers, including autoclave paper pouches, are not acceptable for use with the STERRAD® System because cellulosic materials neutralize the sterilizing agent.

ASP has developed a complete range of sterilization pouches and rolls using Tyvek® 4057B for the STERRAD® System.

More information about the STERRAD® System, including cycle times and performance details, can be found on its **website**.

Other Low-Temperature Oxidative Methods

In addition to hydrogen peroxide gas plasma sterilization, other low-temperature oxidative sterilization methods such as vaporized hydrogen peroxide, nitrogen dioxide and low-temperature peracetic acid vapor gas plasma sterilization have been introduced. DuPont does not have test data for the performance of Tyvek® in each of these low-temperature oxidative sterilization processes; however, the manufacturers of these sterilizers usually have data on Tyvek® compatibility. Please contact the sterilizer manufacturer for more information.

Steam

Tyvek® has been shown to meet packaging criteria for **steam sterilization** under controlled conditions (250°F to 260°F [121°C to 127°C] at 30 psi for 30 minutes). Tyvek® continues to be superior to medical-grade paper when strong, low-linting packaging is required. Tyvek® retains its dimensional stability and integrity with no discoloration when steam sterilized under the controlled conditions mentioned above. Rigid or semi-rigid trays restrict potential shrinkage and wrinkling, which can result in a smoother/tighter lid. Shrinkage of Tyvek® after steam sterilization is less than 1.6%. Tensile strength, microbial barrier and Gurley Hill porosity of Tyvek® are maintained after steam sterilization under the controlled conditions stated above.

Biocompatibility

Biological evaluation of Tyvek® styles for medical and pharmaceutical packaging was performed using testing methodologies according to ISO 10993 and United States Pharmacopeia (USP). In all cases, the styles met all the acceptable performance criteria. This also held true for samples of Tyvek® tested after exposure to sterilization by the EO, gamma and electron-beam processes, proving that Tyvek® meets all the acceptable performance criteria after sterilization. The results of the testing indicate biocompatibility—even after sterilization.

TOXICOLOGICAL RESULTS FOR DUPONT™ TYVEK®

TEST PERFORMED	UNEXPOSED	ETHYLENE OXIDE (EO)	GAMMA IRRADIATION (25 KGY & 50 KGY)	ELECTRON-BEAM (25 KGY & 50 KGY)	STERRAD®
Determination of Extractives from Olefin Polymers ¹		Below maximum allowable percentage			
Hemolysis-Rabbit Blood-ISO ^{2,3}		Non-hemolytic			
L929 MEM Elution Test-USP ⁴		Non-cytotoxic			
ISO-Rabbit Pyrogen Test (Material Mediated) ^{5,6}		Non-pyrogenic			
Kligman Maximization Test-ISO (CSO and NaCl extracts) ^{7,8}		Non-allergenic			
Systemic Injection Test-ISO ⁵		No biological reaction			
Primary Skin Irritation Test-ISO ⁹		Non-irritant			
Short Term Intramuscular Implantation Test-ISO (14 and 28 days) ^{7,8}		Non-irritant			
USP Class VI Test ⁹		0% sensitization			

Tests were based on the following references:
1. 21 CFR 177.1520, Olefin Polymers, Federal Register, Title 21, Chapter 1, 1997.
2. Biological Evaluation of Medical Devices-Part 12: Sample Preparation and Reference Materials, EN/ISO 10993-12, 1997.
3. Biological Evaluation of Medical Devices-Part 4: Selection of Tests for Interactions with Blood, ISO 10993-4, 1992.
4. United States Pharmacopeia 25, National Formulary 20, 2002, <87> Biological Reactivity Tests, In Vitro.
5. Biological Evaluation of Medical Devices-Part 11: Tests for Systemic Toxicity, EN/ISO 10993-1995.
6. Biological Evaluation of Medical Devices-Part 10: Tests for Irritation and Sensitization, EN/ISO 10993-10, 1996.
7. Biological Evaluation of Medical Devices-Part 6: Tests for Local Effects After Implantation, ISO 10993-6, 1995.
8. ASTM Section 13, Volume 13.01 Medical Devices, Designation: F 981-93, 1994.
9. United States Pharmacopeia 25, National Formulary 20, 2002, <88> Biological Reactivity Tests, In Vivo

Featured brands	Our Company	Our Solutions	Support	Ethics & Compliance
Kevlar®	About	Product Lines	Help Center	Position Statements
Nomex®	Careers	Industries	SDS Finder	Ethics Hotline
Corian®	Newsroom	Brands	Brand Licensing	Code of Conduct
Tyvek®	Sustainability	Solution Finder	Supplier Center	REACH
Sorona®	Global Locations			
Danisco®	Investors			
GREAT STUFF™				
Styrofoam™				

